



Advanced
Endoscopy
Devices

NAV SURG ENDOCAM[®]

Navigator Surgical Endoscopy Camera



NSE3500 Micro Endoscopic Video System

- HD Imaging with High Resolution Images
- Light Weight, Flexible, Wireless
- Articulating Display, Easy to View
- Large Storage Capacity
- Eco-friendly Rechargeable Battery
- Excellent for All Endoscopy Procedures
- Easily Adapts to Endoscopes

Optical Zoom	C-Mount	Exposure	AE
Display Screen	3.5", 480*RGB*800	Storage Medium	Max.32GB
Tilt Range	0°~120°	External Power	AC100~240v 50/60Hz DC5v
Image Sensor	2 million pixel high performance sensor	Battery Capacity	1,900mAh 3.7v (Rechargeable & removable, Li-ion polymer battery)
Image Format	JPG	Operating time	Exceed 4hrs
Video Format	Mov, 1080P		
White Balance	AWB		

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md