

Max LAP[®] Full Screen 5mm Laparoscope


ENDOSCOPE FEATURES:

- Bright Sharp Images throughout the Field-Of-View
- Wide Angle Lenses for Field-Of-View of 85 Degrees
- Less Traumatic than a 10mm Scope
- Full Screen Image - No Need to Zoom Video Camera to Fill the Screen
- Ease and Convenience of Use as a Result of High Tech, Lightweight Construction
- Tri-Opt[®] Precision Rod Lens System Manufactured in Germany
- State of the Art Compact Objective System
- High Aperture Fiberoptics Light Transmission
- Multiple Adaptors on Lightpost for use with Various Fiberoptic Light Cables
- Color Coded for Quick Identification
- Ease of Reprocessing: Cold Soak - ETO Gas - Steam Autoclavable
- Endoscopes are warranted to be free from defects in material and workmanship for 3 years from date of purchase

5715B - 5.2mm x 30 degree
Full Screen Wide
Color Code: **Red**

5715F - 5.2mm x 45 degree
Full Screen Wide
Color Code: **Black**

5715A - 5.2mm x 0 degree
Full Screen Wide
Color Code: **Green**



**Advanced
Endoscopy
Devices**

Endoscopy Specialists since 1985

22134 Sherman Way
Canoga Park, CA 91303 - U.S.A.
phone: 818-227-2720 - fax: 818-227-2724
www.aed.md
E-mail: advendodev@earthlink.net