



Advanced  
Endoscopy  
Devices

# Solar 2000

## *Portable LED Lightsource*



### Model # Solar 2000

- Length 12cm
- 4 Hours of Battery Life
- Charging Time: 2 Hours
- Power: 5 watts
- Luminance: 250000 Lux
- Includes LI-ION Charger and Battery

**1-818-227-2720**

Advanced Endoscopy Devices, Inc.  
22134 Sherman Way, Canoga Park, CA 91303 USA  
[www.aed.md](http://www.aed.md) [sales@aed.md](mailto:sales@aed.md)